

## COLLECTION OF FIREBLIGHT SYMPTOMATIC TISSUE FROM APPLE TREES

- If available, use bacterial ooze directly for sample preparation (Picture 1).
- From young twigs, collect 1" sample tissue, taking both healthy and diseased tissue from the transition zone (Picture 2), and slice the shoot longitudinally in half.
- For older bark and active cankers, skim off bark using a knife, need only about 1" length exposed, leaving the cambium layer (Picture 3).

### COLLECTION OF SAMPLE

- If cambium layer looks dry, add a drop of sterile water onto exposed cambium area.
- Rub the exposed cambium layer with a clean disposable pipet tip 3 times or gently touch pipette tip to the ooze, dip pipette into a 200µl microfuge tube containing 25µl sterile water, repeat the process 3 times.

**Note:** a discolored sample gives a poor result.

### TURNING ON AND CONNECTING THE BIORANGER AND TABLET

- Fully charge tablet and BioRanger.
- Turn on tablet and open BioRanger App.
- Turn on BioRanger device and **wait for green blinking light** on top of device.
- Connect tablet and BioRanger device using Connect menu on the tablet (blue light)
- Tablet top right will display 'connected'.
- Select Analysis or Options/Edit methods.
- Analysis method will be indicated in top center.

## BACTERIAL CELL LYSIS

- Collected sample needs to be 'lysed' by heating for 5 minutes at 95°C in the BioRanger.
- Place the 25 µl sample (up to 8 tubes at a time) in the BioRanger.
- Select or Create Lysis Analysis method.
- Options/Edit methods/Add new method/Confirm/Exit.
- Check Lysis step, Set to 95°C for 5 minutes.
- Select Start method, upon start the light on the BioRanger will turn red.
- Status will display in upper left corner of tablet.
- Check lanes you desire to see using the BioRanger App.
- When reaction completed select Options/Quit/Yes.
- When completed cool samples in ice for 5 minutes and let BioRanger rest 5 minutes between runs.

### MIXING SAMPLE WITH OTHER REAGENTS

**Note:** Keep solutions on Ice.

- In new clean tubes, mix the following for each sample needed :  
15µl Optigene mix ISO-001 (Pro-Lab Diagnostics Inc., Round Rock, TX), 2µl DS Primers, 3µl sterile water, 5µl lysed sample.
- Mix four components together by thumping with fingers, let them settle or spin with makeshift centrifuge (Picture 4).
- Run a positive control with pathogen DNA.
- Run a negative control (water, or uninfected plant material).

## RUN ANALYSIS PROTOCOL

- Turn on/connect BioRanger and Tablet, and select analysis method or Options/Edit/Methods Settings: (check boxes) in the BioRanger App.
- Setup, isothermal reaction at 65°C for 15-30 minutes, denaturation at 80°C for 5 minutes, and fluorescence read interval every 30 seconds.
- Other options are, reaction mechanism as LAMP, and Classification approach as Auto
- Afterwards, place tubes with prepared sample/reaction mix in BioRanger and select start method. Indicator light will turn red.
- Make sure SAVE Data is checked.
- Label the lanes with sample information and 'Accept Labels'.
- Select positive and negative labels as preinstalled.
- BioRanger will beep when done and status will display in top left corner of tablet.

### INTERPRETING THE RESULTS

- Watch the results in the BioRanger app on the tablet for real time results.
- A sigmoid curve (S) of positive samples should be observed within 5-10 minutes of the run showing increase in the fluorescence as a peak.

**Note:**

- Letter colors of the sample refer to same color lines on the graph.
- Irregular peaks or delayed peaks can also be seen in in negative samples. Therefore, look at your positive and negative controls and compare them to your sample results.

**Supplies needed:** Clean bench space, ethanol 70%, wipes, ice/cooler, 200ul tubes, pipettes, knife, ISO-001 mix ([www.pro-lab-direct.com](http://www.pro-lab-direct.com)), primer mix, sample, sterile water, BioRanger and an android tablet with BioRanger app installed ([diagenetix.com](http://diagenetix.com)).



Picture 1. Bacterial ooze.




Picture 2. Transition zone between healthy and infected tissue.



Picture 3. Scraped bark to expose cambium layer of fire blight canker.



Picture 4: Additional supplies needed to run loop mediated isothermal DNA amplification



Picture 5: Sigmoid curves for test samples diluted at different concentrations are in gray, violet, blue, yellow, green and blue. Whereas, red line represents the negative sample.

**Note:**

- Size of Threshold peak value (default 10000RFU) will be dependent on the run.
- F=fluorescence, T=temperature, t=time.
- $10^2$  is the detection limit for this kit.

# Detection of Fire Blight *Erwinia amylovora* Pathogen with LOOP MEDIATED ISOTHERMAL DNA AMPLIFICATION

USING THE BIORANGER® BY DIAGENETIX, INC.



Picture 6: BioRanger device (top) and an android Tablet (bottom)

ADDRESS

Dr. Awais Khan  
112 Barton Lab  
15 Castle Creek Drive  
Geneva, NY 14456

CONTACT US

Cornell AgriTech | 315-7872446 |  
[awais.khan@cornell.edu](mailto:awais.khan@cornell.edu)  
<http://blogs.cornell.edu/Khanlab>

